

Kif1a Antibody



PACO27797

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Motor for anterograde axonal transport of synaptic vesicle precursors.

Gene ID:

Kif1a

Uniprot

P33173

Synonyms:

Kinesin-like protein KIF1A (Axonal transporter of synaptic vesicles), Kif1a, Atsv Kif1

Immunogen:

Recombinant Mouse Kinesin-like protein KIF1A protein (1-361AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A