

Kif1c Antibody



PACO27801

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable motor protein.

Gene ID:

Kif1c

Uniprot

O35787

Synonyms:

Kinesin-like protein KIF1C (Kinesin-like protein KIF1D), Kif1c, Kif1d

Immunogen:

Recombinant Rat Kinesin-like protein KIF1C protein (1-272AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A