

Kif1c Antibody, Biotin conjugated



PACO27804

Product Information

Size:	Protein Background:
50ug	Probable motor protein.
Reactivity:	Gene ID:
Rat	Kif1c
Source:	Uniprot
Rabbit	O35787
Isotype:	Synonyms:
IgG	Kinesin-like protein KIF1C (Kinesin-like protein KIF1D), Kif1c, Kif1d
Applications:	Immunogen:
ELISA	Recombinant Rat Kinesin-like protein KIF1C protein (1-272AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A