

KLC4 Antibody



PACO27805

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Kinesin is a microtubule-associated force-producing protein that may play a role in organelle transport. The light chain may function in coupling of cargo to the heavy chain or in the modulation of its ATPase activity.

Gene ID:

KLC4

Uniprot

Q9NSK0

Synonyms:

Kinesin light chain 4 (KLC 4) (Kinesin-like protein 8), KLC4, KNSL8

Immunogen:

Recombinant Human Kinesin light chain 4 protein (401-619AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A