

KRT10 Antibody, FITC conjugated



PACO27847

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KRT10

Uniprot

P13645

Synonyms:

Keratin, type I cytoskeletal 10 (Cytokeratin-10) (CK-10) (Keratin-10) (K10), KRT10, KPP

Immunogen:

Recombinant Human Keratin, type I cytoskeletal 10 protein (1-584AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A