

KRT17 Antibody, FITC conjugated



PACO27855

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Bovine	KRT17
Source:	Uniprot
Rabbit	A1L595
Isotype:	Synonyms:
IgG	Keratin, type I cytoskeletal 17 (Cytokeratin-17) (CK-17) (Keratin-17) (K17), KRT17
Applications:	Immunogen:
ELISA	Recombinant Bovine Keratin, type I cytoskeletal 17 protein (1-441AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A