

L2HGDH Antibody, Biotin conjugated



PACO27888

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

L2HGDH

Uniprot

A7MBI3

Synonyms:

L-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99.2) (Duranin), L2HGDH

Immunogen:

Recombinant Bovine L-2-hydroxyglutarate dehydrogenase, mitochondrial protein (53-463AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A