

LHFPL6 Antibody



PACO27949

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LHFPL6

Uniprot

Q9Y693

Synonyms:

LHFPL tetraspan subfamily member 6 protein (Lipoma HMGIC fusion partner), LHFPL6, LHFPL

Immunogen:

Recombinant Human LHFPL tetraspan subfamily member 6 protein (31-83AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A