

LILRA4 Antibody, Biotin conjugated



PACO27964

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as receptor for class I MHC antigens. Ligand binding leads to the activation of eosinophils and to the release of RNASE2, IL4 and leukotriene C4.

Gene ID:

LILRA4

Uniprot

P59901

Synonyms:

Leukocyte immunoglobulin-like receptor subfamily A member 4 (CD85 antigen-like family member G) (Immunoglobulin-like transcript 7) (ILT-7) (CD antigen CD85g), LILRA4, ILT7

Immunogen:

Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 4 protein (289-424AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A