

LRRN2 Antibody



PACO28025

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

receptor activity, cell adhesion, signal transduction.

Gene ID:

LRRN2

Uniprot

O75325

Synonyms:

Leucine-rich repeat neuronal protein 2 (Glioma amplified on chromosome 1 protein)
(Leucine-rich repeat neuronal protein 5), LRRN2, GAC1 LRRN5


Immunogen:

Recombinant Human Leucine-rich repeat neuronal protein 2 protein (201-500AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO28025 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO28025 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).