

LSM1 Antibody, HRP conjugated



PACO28034

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in replication-dependent histone mRNA degradation. Binds specifically to the 3'-terminal U-tract of U6 snRNA.

Gene ID:

LSM1

Uniprot

O15116

Synonyms:

U6 snRNA-associated Sm-like protein LSm1 (Cancer-associated Sm-like) (Small nuclear ribonuclear CaSm), LSM1, CASM

Immunogen:

Recombinant Human U6 snRNA-associated Sm-like protein LSm1 protein (1-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A