

# LSM3 Antibody, FITC conjugated



PACO28043

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds specifically to the 3'-terminal U-tract of U6 snRNA.

**Gene ID:**

LSM3

**Uniprot**

P62310

**Synonyms:**

U6 snRNA-associated Sm-like protein LSm3, LSM3

**Immunogen:**

Recombinant Human U6 snRNA-associated Sm-like protein LSm3 protein (2-102AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A