

LSM4 Antibody, HRP conjugated



PACO28046

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds specifically to the 3'-terminal U-tract of U6 snRNA.

Gene ID:

LSM4

Uniprot

Q9Y4Z0

Synonyms:

U6 snRNA-associated Sm-like protein LSm4 (Glycine-rich protein) (GRP), LSM4

Immunogen:

Recombinant Human U6 snRNA-associated Sm-like protein LSm4 protein (1-139AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A