

# LY6G6C Antibody, HRP conjugated



PACO28059

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LY6G6C

**Uniprot**

O95867

**Synonyms:**

Lymphocyte antigen 6 complex locus protein G6c, LY6G6C, C6orf24 G6C NG24

**Immunogen:**

Recombinant Human Lymphocyte antigen 6 complex locus protein G6c protein (19-99AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A