

LY6G6C Antibody, Biotin conjugated



PACO28061

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LY6G6C

Uniprot

O95867

Synonyms:

Lymphocyte antigen 6 complex locus protein G6c, LY6G6C, C6orf24 G6C NG24

Immunogen:

Recombinant Human Lymphocyte antigen 6 complex locus protein G6c protein (19-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A