

Ly75 Antibody, FITC conjugated



PACO28072

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as an endocytic receptor to direct captured antigens from the extracellular space to a specialized antigen-processing compartment. Causes reduced proliferation of B lymphocytes.

Gene ID:

Ly75

Uniprot

Q60767

Synonyms:

Lymphocyte antigen 75 (Ly-75) (DEC-205) (CD antigen CD205), Ly75, Cd205

Immunogen:

Recombinant Mouse Lymphocyte antigen 75 protein (28-211AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A