

# Ly75 Antibody, Biotin conjugated



PACO28073

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Acts as an endocytic receptor to direct captured antigens from the extracellular space to a specialized antigen-processing compartment. Causes reduced proliferation of B lymphocytes.

**Gene ID:**

Ly75

**Uniprot**

Q60767

**Synonyms:**

Lymphocyte antigen 75 (Ly-75) (DEC-205) (CD antigen CD205), Ly75, Cd205

**Immunogen:**

Recombinant Mouse Lymphocyte antigen 75 protein (28-211AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A