

LY96 Antibody, HRP conjugated



PACO28075

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds bacterial lipopolysaccharide (LPS). Cooperates with TLR4 in the innate immune response to bacterial lipopolysaccharide (LPS), and with TLR2 in the response to cell wall components from Gram-positive and Gram-negative bacteria. Enhances TLR4-dependent activation of NF-kappa-B. Cells expressing both LY96 and TLR4, but not TLR4 alone, respond to LPS.

Gene ID:

LY96

Uniprot

Q9Y6Y9

Synonyms:

Lymphocyte antigen 96 (Ly-96) (ESOP-1) (Protein MD-2), LY96, ESOP1 MD2

Immunogen:

Recombinant Human Lymphocyte antigen 96 protein (19-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A