

SNX3 Antibody, Biotin conjugated



PACO28081

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Phosphoinositide-binding protein required for multivesicular body formation. Specifically binds phosphatidylinositol-3-phosphate (PtdIns(P3)). Plays a role in protein transport between cellular compartments. Promotes stability and cell surface expression of epithelial sodium channel (ENAC) subunits SCNN1A and SCNN1G By similarity. Not involved in EGFR degradation.

Gene ID:

SNX3

Uniprot

O60493

Synonyms:

Sorting nexin-3 (Protein SDP3), SNX3

Immunogen:

Recombinant Human Sorting nexin-3 protein (3-162AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A