

LYRM1 Antibody



PACO28094

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May promote cell proliferation and inhibition of apoptosis of preadipocytes.

Gene ID:

LYRM1

Uniprot

O43325

Synonyms:

LYR motif-containing protein 1, LYRM1


Immunogen:

Recombinant Human LYR motif-containing protein 1 protein (1-122AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO28094 at dilution of 1:100.