

# LYRM1 Antibody, Biotin conjugated



PACO28097

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May promote cell proliferation and inhibition of apoptosis of preadipocytes.

**Gene ID:**

LYRM1

**Uniprot**

O43325

**Synonyms:**

LYR motif-containing protein 1, LYRM1

**Immunogen:**

Recombinant Human LYR motif-containing protein 1 protein (1-122AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A