

# MAGEA10 Antibody, FITC conjugated



PACO28104

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

**Gene ID:**

MAGEA10

**Uniprot**

P43363

**Synonyms:**

Melanoma-associated antigen 10 (Cancer/testis antigen 1.10) (CT1.10) (MAGE-10 antigen), MAGEA10, MAGE10

**Immunogen:**

Recombinant Human Melanoma-associated antigen 10 protein (2-225AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A