

MAGEA10 Antibody, Biotin conjugated



PACO28105

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

Gene ID:

MAGEA10

Uniprot

P43363

Synonyms:

Melanoma-associated antigen 10 (Cancer/testis antigen 1.10) (CT1.10) (MAGE-10 antigen), MAGEA10, MAGE10

Immunogen:

Recombinant Human Melanoma-associated antigen 10 protein (2-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A