

MAGEA8 Antibody, Biotin conjugated



PACO28117

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

Gene ID:

MAGEA8

Uniprot

P43361

Synonyms:

Melanoma-associated antigen 8 (Cancer/testis antigen 1.8) (CT1.8) (MAGE-8 antigen), MAGEA8, MAGE8

Immunogen:

Recombinant Human Melanoma-associated antigen 8 protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A