

MAGEA9 Antibody



PACO28118

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

Gene ID:

MAGEA9

Uniprot

P43362

Synonyms:

Melanoma-associated antigen 9 (Cancer/testis antigen 1.9) (CT1.9) (MAGE-9 antigen), MAGEA9; MAGEA9B, MAGE9 MAGEA9A;


Immunogen:

Recombinant Human Melanoma-associated antigen 9 protein (1-315AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO28118 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung cancer using PACO28118 at dilution of 1:100.