

PSMA7 Antibody, HRP conjugated



PACO28123

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PSMA7

Uniprot

O14818

Synonyms:

Proteasome subunit alpha type-7 (EC 3.4.25.1) (Proteasome subunit RC6-1) (Proteasome subunit XAPC7), PSMA7, HSPC

Immunogen:

Recombinant Human Proteasome subunit α ; type-7 protein (5-232AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A