

MFAP3 Antibody



PACO28214

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the elastin-associated microfibrils.

Gene ID:

MFAP3

Uniprot

P55082

Synonyms:

Microfibril-associated glycoprotein 3, MFAP3

Immunogen:

Recombinant Human Microfibril-associated glycoprotein 3 protein (171-362AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A