

MFAP3L Antibody



PACO28218

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

cytoplasm, nucleus, plasma membrane.

Gene ID:

MFAP3L

Uniprot

O75121

Synonyms:

Microfibrillar-associated protein 3-like (Testis development protein NYD-SP9), MFAP3L, KIAA0626


Immunogen:

Recombinant Human Microfibrillar-associated protein 3-like protein (29-102AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO28218 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO28218 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human testis tissue using PACO28218 at dilution of 1:100.