

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions in replication-dependent translation of histone mRNAs which differ from other eukaryotic mRNAs in that they do not end with a poly-A tail but a stem-loop. May participate in circularizing those mRNAs specifically enhancing their translation.

Gene ID:

MIF4GD

Uniprot

A9UHW6

Synonyms:

MIF4G domain-containing protein (SLBP-interacting protein 1) (hSLIP1), MIF4GD, SLIP1

Immunogen:

Recombinant Human MIF4G domain-containing protein (1-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A