

MPPED1 Antibody, Biotin conjugated



PACO28317

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May have metallophosphoesterase activity (in vitro).

Gene ID:

MPPED1

Uniprot

O15442

Synonyms:

Metallophosphoesterase domain-containing protein 1 (EC 3.1. -. -) (Adult brain protein 239) (239AB), MPPED1, C22orf1 FAM1A

Immunogen:

Recombinant Human Metallophosphoesterase domain-containing protein 1 protein (1-326AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A