

MPZL2 Antibody, FITC conjugated



PACO28324

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Mediates homophilic cell-cell adhesion.

Gene ID:

MPZL2

Uniprot

O60487

Synonyms:

Myelin protein zero-like protein 2 (Epithelial V-like antigen 1), MPZL2, EVA EVA1

Immunogen:

Recombinant Human Myelin protein zero-like protein 2 protein (30-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A