

MRPL42 Antibody, Biotin conjugated



PACO28357

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MRPL42

Uniprot

Q9Y6G3

Synonyms:

39S ribosomal protein L42, mitochondrial (L42mt) (MRP-L42) (39S ribosomal protein L31, mitochondrial) (L31mt) (MRP-L31) (Mitochondrial large ribosomal subunit protein mL42), MRPL42, MRPL31 MRPS32 RPML31

Immunogen:

Recombinant Human 39S ribosomal protein L42, mitochondrial protein (33-142AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A