

MRPL54 Antibody, FITC conjugated



PACO28364

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MRPL54

Uniprot

Q6P161

Synonyms:

39S ribosomal protein L54, mitochondrial (L54mt) (MRP-L54) (Mitochondrial large ribosomal subunit protein mL54), MRPL54

Immunogen:

Recombinant Human 39S ribosomal protein L54, mitochondrial protein (15-138AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A