

# MRPL55 Antibody, FITC conjugated



PACO28368

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of the mitochondrial ribosome large subunit (39S) which comprises a 16S rRNA and about 50 distinct proteins.

**Gene ID:**

MRPL55

**Uniprot**

Q7Z7F7

**Synonyms:**

39S ribosomal protein L55, mitochondrial (L55mt) (MRP-L55) (Mitochondrial large ribosomal subunit protein bL31m) (Mitochondrial large ribosomal subunit protein mL55), MRPL55

**Immunogen:**

Recombinant Human 39S ribosomal protein L55, mitochondrial protein (34-128AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A