

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

mitochondrial inner membrane, mitochondrial small ribosomal subunit, structural constituent of ribosome, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization, translation.

**Gene ID:**

MRPS18C

**Uniprot**

Q9Y3D5

**Synonyms:**

28S ribosomal protein S18c, mitochondrial (MRP-S18-c) (Mrps18-c) (S18mt-c) (28S ribosomal protein S18-1, mitochondrial) (MRP-S18-1) (Mitochondrial small ribosomal subunit protein bS18c) (Mitochondrial small ribosomal subunit protein bS18m), MRPS18C

**Immunogen:**

Recombinant Human 28S ribosomal protein S18c, mitochondrial protein (1-142AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A