

Ms4a1 Antibody, Biotin conjugated



PACO28397

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This protein may be involved in the regulation of B-cell activation and proliferation.

Gene ID:

Ms4a1

Uniprot

P19437

Synonyms:

B-lymphocyte antigen CD20 (B-cell differentiation antigen Ly-44) (Lymphocyte antigen 44) (Membrane-spanning 4-domains subfamily A member 1) (CD antigen CD20), Ms4a1, Cd20 Ly-44 Ms4a2

Immunogen:

Recombinant Mouse B-lymphocyte antigen CD20 protein (111-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A