

# Msln Antibody, HRP conjugated



PACO28407

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Membrane-anchored forms may play a role in cellular adhesion. Megakaryocyte-potentiating factor (MPF) may potentiate megakaryocyte colony formation.

**Gene ID:**

Msln

**Uniprot**

Q61468

**Synonyms:**

Mesothelin (Pre-pro-megakaryocyte-potentiating factor) [Cleaved into: Megakaryocyte-potentiating factor (MPF); Mesothelin, cleaved form], Msln, Mes Mpf

**Immunogen:**

Recombinant Mouse Mesothelin protein (298-600AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A