

# MUC2 Antibody, HRP conjugated



PACO28431

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Coats the epithelia of the intestines, airways, and other mucus membrane-containing organs. Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces. Major constituent of both the inner and outer mucus layers of the colon and may play a role in excluding bacteria from the inner mucus layer.

**Gene ID:**

MUC2

**Uniprot**

Q02817

**Synonyms:**

Mucin-2 (MUC-2) (Intestinal mucin-2), MUC2, SMUC

**Immunogen:**

Recombinant Human Mucin-2 protein (36-240AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A