

Myc Antibody, HRP conjugated



PACO28455

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Myc

Uniprot

P01108

Synonyms:

Myc proto-oncogene protein (Proto-oncogene c-Myc) (Transcription factor p64), Myc

Immunogen:

Recombinant Mouse Myc proto-oncogene protein (1-439AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A