

PACO28470

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

nucleus, transcription factor activity, sequence-specific DNA binding, zinc ion binding, regulation of transcription, DNA-templated, transcription, DNA-templated.

Gene ID:

MYNN

Uniprot

Q9NPC7

Synonyms:

Myoneurin (Zinc finger and BTB domain-containing protein 31), MYNN, OSZF ZBTB31


Immunogen:

Recombinant Human Myoneurin protein (51-330AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO28470 at dilution of 1:100.