

# MYNN Antibody, Biotin conjugated



PACO28473

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

nucleus, transcription factor activity, sequence-specific DNA binding, zinc ion binding, regulation of transcription, DNA-templated, transcription, DNA-templated.

**Gene ID:**

MYNN

**Uniprot**

Q9NPC7

**Synonyms:**

Myoneurin (Zinc finger and BTB domain-containing protein 31), MYNN, OSZF ZBTB31

**Immunogen:**

Recombinant Human Myoneurin protein (51-330AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A