

# NADK Antibody, Biotin conjugated



PACO28481

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytosol, NAD+ kinase activity, ATP metabolic process, NAD metabolic process, phosphorylation, small molecule metabolic process, vitamin metabolic process, water-soluble vitamin metabolic process.

**Gene ID:**

NADK

**Uniprot**

O95544

**Synonyms:**

NAD kinase (EC 2.7.1.23) (Poly(P))/ATP NAD kinase, NADK

**Immunogen:**

Recombinant Human NAD kinase protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A