

NAPA Antibody, Biotin conjugated



PACO28489

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for vesicular transport between the endoplasmic reticulum and the Golgi apparatus.

Gene ID:

NAPA

Uniprot

P54920

Synonyms:

Alpha-soluble NSF attachment protein (SNAP-alpha) (N-ethylmaleimide-sensitive factor attachment protein alpha), NAPA, SNAPA

Immunogen:

Recombinant Human Alpha-soluble NSF attachment protein (1-295AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A