

NAT2 Antibody, HRP conjugated



PACO28503

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Participates in the detoxification of a plethora of hydrazine and arylamine drugs. Catalyzes the N- or O-acetylation of various arylamine and heterocyclic amine substrates and is able to bioactivate several known carcinogens.

Gene ID:

NAT2

Uniprot

P11245

Synonyms:

Arylamine N-acetyltransferase 2 (EC 2.3.1.5) (Arylamide acetylase 2) (N-acetyltransferase type 2) (NAT-2) (Polymorphic arylamine N-acetyltransferase) (PNAT), NAT2, AAC2

Immunogen:

Recombinant Human Arylamine N-acetyltransferase 2 protein (1-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A