

NDUFA3 Antibody, Biotin conjugated



PACO28521

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NDUFA3

Uniprot

O95167

Synonyms:

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 (Complex I-B9) (CI-B9) (NADH-ubiquinone oxidoreductase B9 subunit), NDUFA3

Immunogen:

Recombinant Human NADH dehydrogenase [ubiquinone] 1 α subcomplex subunit 3 protein (2-84AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A