

# NFE2L1 Antibody, HRP conjugated



PACO28551

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Activates erythroid-specific, globin gene expression.

**Gene ID:**

NFE2L1

**Uniprot**

Q14494

**Synonyms:**

Endoplasmic reticulum membrane sensor NFE2L1 (Locus control region-factor 1) (LCR-F1) (Nuclear factor erythroid 2-related factor 1) (NF-E2-related factor 1) (NFE2-related factor 1) (Nuclear factor, erythroid derived 2, like 1) (Protein NRF1, p120 form) (Transcription factor 11) (TCF-11) [Cleaved into: Transcription factor NRF1 (Protein NRF1, p110 form)], NFE2L1, HBZ17 NRF1 TCF11

**Immunogen:**

Recombinant Human Endoplasmic reticulum membrane sensor NFE2L1 protein (251-550AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A