

NFE2L1 Antibody, Biotin conjugated



PACO28553

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Activates erythroid-specific, globin gene expression.

Gene ID:

NFE2L1

Uniprot

Q14494

Synonyms:

Endoplasmic reticulum membrane sensor NFE2L1 (Locus control region-factor 1) (LCR-F1) (Nuclear factor erythroid 2-related factor 1) (NF-E2-related factor 1) (NFE2-related factor 1) (Nuclear factor, erythroid derived 2, like 1) (Protein NRF1, p120 form) (Transcription factor 11) (TCF-11) [Cleaved into: Transcription factor NRF1 (Protein NRF1, p110 form)], NFE2L1, HBZ17 NRF1 TCF11

Immunogen:

Recombinant Human Endoplasmic reticulum membrane sensor NFE2L1 protein (251-550AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A