

Nid1 Antibody, FITC conjugated



PACO28568

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sulfated glycoprotein widely distributed in basement membranes and tightly associated with laminin. Also binds to collagen IV and perlecan. It probably has a role in cell-extracellular matrix interactions.

Gene ID:

Nid1

Uniprot

P10493

Synonyms:

Nidogen-1 (NID-1) (Entactin), Nid1, Ent

Immunogen:

Recombinant Mouse Nidogen-1 protein (428-665AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A