

# Nid1 Antibody, Biotin conjugated



PACO28569

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Sulfated glycoprotein widely distributed in basement membranes and tightly associated with laminin. Also binds to collagen IV and perlecan. It probably has a role in cell-extracellular matrix interactions.

**Gene ID:**

Nid1

**Uniprot**

P10493

**Synonyms:**

Nidogen-1 (NID-1) (Entactin), Nid1, Ent

**Immunogen:**

Recombinant Mouse Nidogen-1 protein (428-665AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A