

# Nkx2-2 Antibody, HRP conjugated



PACO28575

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Acts as a transcriptional activator. Required for the maintenance of NEUROD1 expression in the hormone-producing endocrine cells of the pancreas. May be involved in specifying diencephalic neuromeric boundaries, and in controlling the expression of genes that play a role in axonal guidance. Associates with chromatin at the NEUROD1 promoter region. Binds to a subset of consensus elements within the NEUROD1 promoter.

**Gene ID:**

Nkx2-2

**Uniprot**

P42586

**Synonyms:**

Homeobox protein Nkx-2.2 (Homeobox protein NK-2 homolog B), Nkx2-2, Nkx-2.2 Nkx2b

**Immunogen:**

Recombinant Mouse Homeobox protein Nkx-2.2 protein (1-273AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A